

JUNO, MIRACURE Infrared Heating Lamp

How to use

1. Turn on power switch
2. One to Two treatment a day, 15 ~ 20 min every time.
3. The distance between the redheating lamp and body , appx 12" ~ 16" .
Note: The distance varies from person to person, Adjust accordingly to fit desired comfort level for treatment.
4. Ideal temperature when redheating lamp hits skin surface 115° ~ 130°F.
5. Long term use increase better heath benefits.


How to Remove the FOAM



Put cross screwdriver in the hole, and take out 4 screws, and Remove the foam.

- The Infrared Heating Lamp is effective when deep heating of body tissue is desired.
- Beneficial effects include decrease joint stiffness, increased vasodilation, muscle spasm relief and reduced pain from ligamentous sprains and strains.
- It works by increasing the local circulation of blood, lymph and the activity of the nervous system.
- Decrease painful edema
- Reduces pain in joints, muscles, bones
- Relaxes muscle armoring and spasms
- Improves microcirculation

JUNO, MIRACURE Infrared Heating Lamp



Caution

1. Person suffering with eye complications should not look directly into the redheating lamp, Both eyes should be covered with a towel.
2. Avoid exposing any open wounds or cuts to the redheating lamp.
3. Person(s) pregnant should consult with physician.
4. Don' t use over maximum 30 min.
5. If you have a kind of allegy, please consult with physician.
6. If you feel too much hot, please don' t use it.
7. If you have inflammation on the skin, please consult with physician.
8. Before use infrared heating lamp, please consult with physician.

WARRANTY

Neomedic warrants that Infrared Heating Lamp is free of defects in material and workmanship. This warranty shall remain in effect for on year (12 months) from the date of original consumer purchase. If this product fails to function during the one year warranty period due to a defect in material or workmanship, Company or the selling dealer will repair or replace the product without charge within a period of thirty (30) days from the date on which the product is returned to the company or the dealer. All repairs to the product must be performed by a service center authorized by the company.

Company shall not be liable in any event for incidental or consequential damages.

1. A written claim must be made within the warranty period to the company or the selling dealer.
2. The product must be returned to the company or the selling dealer by the owner along with a copy of the purchase receipt



Toll Free: 866-990-1168

www.air1000.com, info@air1000.com